

AUGUST

WEEK 4

Aim: Check whether proteins are expressing in inclusion bodies or not

Protocol:


- Primary inoculum- 10 ml lb + kan(50ug/mL). Incubate for 12hr
- Secondary inoculum- 1050 ul primary inoculum 105 ml lb+ kan (50ug/mL). Incubate till OD reaches 0.6. (preserve 1ml of secondary inoculum as a control)
- Add 1mM IPTG(200ul of 0.5M IPTG) and incubate at 16°C for 12hr.
- Pellet it down. Take some pellet and add loading dye buffer.
- Take some pellet and add lysis buffer.
- Sonicate it and then centrifuge at 4C for 11000 rpm for 45 mins..
- Dissolve both pellet and supernatant in loading buffer separately.
- Perform SDS [Stacking Gel (5%), Resolving Gel (10%)]

SDS Page

Annotations :

- 1:Sample from Supernatant (After sonication)
- 2:Sample from Pellet (After Sonication)
- 3:Sample from Pellet (After Sonication)
- 4:Sample from Supernatant (After sonication)

Result: The protein is coming in both supernatant and pellet at 16°C



Protein Extraction

- 1) Inoculation
- 2) Induction
- 3) Protein Purification

Inoculation and Induction

- Primary inoculum incubated overnight (10mL) (13hrs).
- Secondary inoculum (1L) incubated till OD600 reached 0.6
- Secondary inoculum induced with 1mM IPTG and incubated overnight at 16 C.
- Cells pelleted down at 13,000 RPM for 15 mins and stored at -20 C.

Lysis and Sonication

27/08/2021

Lysis Buffer

Components	Stock concentration	Working Concentration	Volume required for 50 mL solution
Sodium Phosphate Buffer (pH 7)	0.1 M	50 mM	10 mL
PMSF	1 M	1 mM	20 uL
Sodium Chloride	2 M	0.5 M	5 mL
BME		0.05 % of lysis buffer	10 uL
Glycerol		5 % of lysis buffer	1 mL
Triton X-100	100 %	0.5 % of lysis buffer	100 uL
Protease Inhibitor	1 tablet(EDTA free)	1 tablet (EDTA free)	1 tablet (EDTA free)
Autoclaved Milli Q	To make up 50 mL		3.87 mL

Sonication

On time	30 secs
Off time	30 secs
Amplitude	55%
Total cycles	15
Total time	15 mins

Sonication carried out in pulse mode

Protein Purification (Ni-NTA)

31/08/2021

Equilibration Buffer

Component	Stock Concentration	Working Concentration	Volume for 50 mL
Sodium Phosphate Buffer (pH 7)	0.1 M	50 mM	10 mL
PMSF	1 M	1 mM	20 uL
Sodium Chloride	2 M	0.5 M	5 mL
BME		0.05 %	10 uL
Glycerol		5 %	1 mL
Triton X-100	100 %	0.5 %	100 uL
Autoclaved Milli Q	Make upto 50 mL		3.87 mL

Wash Buffer (50 mL)

Component	Stock Concentration	Working Concentration	Volume for 50 mL Wash Buffer
Tris Base (pH 8)	1 M	20 mM	1 mL
Sodium Chloride	5 M	400 mM	4 mL
EDTA	0.5 M	1 mM	100 μ L
Imidazole	1 M	25 mM	1.25 mL
Autoclaved MilliQ	Make upto 50 mL		43.65 mL

Elution Buffer (50 mL)

Component	Stock Concentration	Working Concentration	Volume for 50 mL Elution Buffer
Tris Base (pH 8)	1 M	20 mM	1 mL
Sodium Chloride	5 M	150 mM	1.5 mL
EDTA	0.5 M	1 mM	100 μ L
Imidazole	1 M	250 mM	12.5 mL
Autoclaved MilliQ	Make upto 50 mL		34.9 mL

Ni-NTA Purification

- The flow through (filtrate collected after loading solubilised protein) and wash buffer were collected separately in 50 mL falcon tube.
- The elution flow through was collected in fractions in 20, 2 mL Eppendorf tubes.