

FUNGUS IDENTIFICATION

WEEK 2

PCR OF 1 SAMPLE

Protocol:

Taq buffer	5 uL
MgCl ₂	1 uL
dNTP	1uL
Taq Polymerase	1 uL
10uM FP	1 uL
10uM RP	1 uL
Template(120 ng/uL)	3 uL
Water	8.5 uL
Total	50 uL

Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
1	165.5	1.66	1.24

NO BAND OBSERVED

PCR CLEAN UP OF SAMPLE 1

PROTOCOL :

1. Add 5 volumes Buffer PB to 1 volume of the PCR reaction and mix. If the color of the mixture is orange or violet, add 10 μ l 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn yellow.
2. Place a QIAquick column in a provided 2 ml collection.
- 3.3. To bind DNA, apply the sample to the QIAquick column and centrifuge for 30–60s.
4. Discard flow-through and place the QIAquick column back in the same tube.
5. To wash, add 600 μ l Buffer PE to the QIAquick column & centrifuge for 30–60 s.
6. Discard flow-through and place the QIAquick column back into the same tube.
7. Centrifuge the QIAquick column once more in the provided 2 ml collection tube for 1 min to remove residual wash buffer.
8. Place each QIAquick column in a clean 1.5 ml microcentrifuge tube.
9. To elute DNA, add 25 μ l Buffer EB (10 mM Tris·Cl, pH 8.5) at 70 degree Celsius or water (pH 7.0–8.5) to the center of the QIAquick membrane let the column stand for 1 min in water bath at 70 degree and centrifuge the column for 1 min. For increased DNA concentration, add 25 μ l elution buffer to the center of the QIAquick membrane, let the column stand for 1 min and then centrifuge.

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
1	24.8	1.29	0.89

PCR OF 1 SAMPLE

Protocol:

Master Mix	10 uL
10mM FP	0.5 uL
10mM RP	0.5uL
Template(100-150) ng/uL	1 uL
Water	8.5 uL
Total	20 uL

Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
1	431.6	1.81	2.13

NO BAND OBSERVED

PCR CLEAN UP OF SAMPLE 1

PROTOCOL :

1. Add 5 volumes Buffer PB to 1 volume of the PCR reaction and mix. If the color of the mixture is orange or violet, add 10 μ l 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn yellow.
2. Place a QIAquick column in a provided 2 ml collection.
3. To bind DNA, apply the sample to the QIAquick column and centrifuge for 30–60s.
4. Discard flow-through and place the QIAquick column back in the same tube.
5. To wash, add 600 μ l Buffer PE to the QIAquick column & centrifuge for 30–60 s.
6. Discard flow-through and place the QIAquick column back into the same tube.
7. Centrifuge the QIAquick column once more in the provided 2 ml collection tube for 1 min to remove residual wash buffer.
8. Place each QIAquick column in a clean 1.5 ml microcentrifuge tube.
9. To elute DNA, add 25 μ l Buffer EB (10 mM Tris·Cl, pH 8.5) at 70 degree Celsius or water (pH 7.0–8.5) to the center of the QIAquick membrane let the column stand for 1 min in water bath at 70 degree and centrifuge the column for 1 min. For increased DNA concentration, add 25 μ l elution buffer to the center of the QIAquick membrane, let the column stand for 1 min and then centrifuge.

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
1	15.7	1.58	1.68

PCR OF 6&7 SAMPLE

Protocol:


Master Mix	12.5 uL
10uM FP	0.625 uL
10uM RP	0.625uL
Template(250 ng/uL)	3 → 1.86 uL 4 → 0.86 uL
Water	3 → 9.39 uLuL
Total	25uL

Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
6	375.5	1.82	2.18
7	379.4	1.83	2.08



PCR OF 6&7 SAMPLE

Protocol:

Master Mix	10 uL
10uM FP	0.5 uL
10uM RP	0.5uL
Template(100-150) ng/uL	1 uL
Water	8.5 uL
Total	20 uL

Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

PCR OF 2&4 SAMPLE

Protocol:

Master Mix	12.5 uL
10uM FP	0.625 uL
10uM RP	0.625uL
Template(250 ng/uL)	2 → 2.2uL 4 → 0.97uL
Water	2 → 9.05 uL 4 → 10.28 uL
Total	25 uL


Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

RESULT

S.No.	Concentration (ng/uL)	A260/A280	A260/A2 30
6	366.4	1.83	2.12
4	425.2	1.83	2.12
2	374	1.84	2.12
7	390.3	1.84	2.08

BAND FOR 2ND SAMPLE IS MISSING



PCR OF 2&4 SAMPLE

Protocol:

Master Mix	10 uL
10uM FP	0.5 uL
10uM RP	0.5uL
Template(100-150) ng/uL	1 uL
Water	8.5 uL
Total	20 uL

Reaction:

95°C	5 min	
94°C	30 s	X 35
50°C	30 s	
72°C	1:30 min	
72°C	10 min	
4°C	hold	

RESULT

S.No	Concentration (ng/uL)	A260/A280	A260/A230
2	375.5	1.82	2.18
4	379.4	1.83	2.08

