

# OD-based growth model

---

## Introduction

Scattering of light by a replicating population over time, and hence the changing OD is a measure to model the growth of a fungal species over time. This experiment aims to determine the growth rate and lag time of fungal species identified in our surroundings.

## Materials

- › Fungal samples cultures on slant PDA(Potato dextrose agar)
- › Potato dextrose broth(PDB)
- › Test tubes
- › Phosphate buffered saline(PBS)
- › Tween 20
- › 10mL Falcon tubes
- › Vortex machine
- › Hemocytometer
- › 96-well plate
- › Plate reader

## Procedure

### Preparation of spore suspension

1. Keep the necessary items in the biosafety cabinet and switch on the UV light for sterilization.
2. Switch off the UV light after 15 mins and wipe off the work area with 70% ethanol.
3. Add 5 mL of PBS to the slant test tube, swirl, and pour the PBS into a falcon tube.
4. Add 0.1% (of total volume of PBS added, 5 $\mu$ L, here) Tween 20 to the tube.
5. Vortex suspension for 15 seconds at 2000 rpm for uniform mixing of spores in the PBS.

### Counting of spores using a hemocytometer

6. Clean the hemocytometer with 70% ethanol and kim wipes and place its coverslip covering both the chambers.
7. Load the spore suspension in the two chamber.

8. Place the hemocytometer under microscope and view it under 10X and 20X to clearly visualise the spores.
9. Select the box to be considered according to the size of the spores and count them.
10. Calculate the number of spores in the entire 5ml of spore suspension.

## Setting up the 96-well plate

11. To keep the number of spores constant in all the wells, dilute the spore suspension using PBS.
12. Prepare dilutions to make suspension of  $4 \times 10^4$  spores/mL.
13. 3 wells(triplicate) are to be used for each fungal species in a 96-well plate. Fill the wells with 100 $\mu$ L of PDB.
14. Add 100 $\mu$ L of spore suspension to make the final count of spores in the well as  $2 \times 10^4$  spores/mL.
15. For blank, add 100 $\mu$ L of PDB and 100 $\mu$ L of PBS.

Well1												
	1	2	3	4	5	6	7	8	9	10	11	12
A	1A		2A		3A		Blank					
B	1A		2A		3A		Blank					
C	1A		2A		3A		Blank					
D												
E												
F	1B		2B		3B		Blank					
G	1B		2B		3B		Blank					
H	1B		2B		3B		Blank					

Note: For our experiment , we used 3 fungal species and made 2 separate slant PDA cultures(A and B) for them

1: *Aspergillus niger*(1A and 1B)

2: *Aspergillus versicolor*(2A and 2B)

3: *Rhizopus oryzae*(3A and 3B)

16. Keep the plate in a shaker cum incubator at 28°C.
17. Take the OD measurements at 605nm in a plate reader for a fixed interval of time.. (For our purpose, we took 4 times a day at an interval of 6 hr for 4-5 days.)
18. Record the reading.

# IC50 of an Antifungal

---

## Introduction

IC<sub>50</sub> of an antifungal here would indicate how much of the chemical would be needed to inhibit, in vitro, the given fungal species. (<https://en.wikipedia.org/wiki/IC50>)

Lower the IC<sub>50</sub>, better the efficiency of the antifungal.

## Materials

- › Fungal samples cultures on slant PDA(Potato dextrose agar)
- › Potato dextrose broth(PDB)
- › Test tubes
- › Phosphate buffered saline(PBS)
- › Tween 20
- › 10mL Falcon tubes
- › Nipagin (antifungal)
- › 100% ethanol
- › Vortex machine
- › Hemocytometer
- › 96-well plate
- › Plate reader
- ›

## Procedure

### Preparation of spore suspension

1. Keep the necessary items in the biosafety cabinet and switch on the UV light for sterilization.
2. Switch off the UV light after 15 mins and wipe off the work area with 70% ethanol.
3. Add 5 mL of PBS to the slant test tube, swirl, and pour the PBS into a falcon tube.
4. Add 0.1% (of total volume of PBS added, 5µL, here) Tween 20 to the tube.
5. Vortex suspension for 15 seconds at 2000 rpm for uniform mixing of spores in the PBS.

## Counting of spores using a hemocytometer

6. Clean the hemocytometer with 70% ethanol and kim wipes and place its coverslip covering both the chambers.
7. Load the spore suspension in the two chamber.
8. Place the hemocytometer under microscope and view it under 10X and 20X to clearly visualise the spores.
9. Select the box to be considered according to the size of the spores and count them.
10. Calculate the number of spores in the entire 5ml of spore suspension.

## Use of antifungal-Niapagin

11. To make 0.1g/mL concentration of Nipagin, add 0.1 g of Nipagin to 100% ethanol.
12. Dilutions can be made according to the need using 100% ethanol.

## Setting up the 96-well plate

13. To keep the number of spores constant in all the wells, dilute the spore suspension using PBS.
14. Prepare dilutions to make suspension of  $4 \times 10^4$  spores/mL and  $8 \times 10^4$  spores/mL.
15. 3 wells(triplicate) are to be used for each concentration of Nipagin in a 96-well plate. Fill the wells with 100 $\mu$ L of PDB.
16. Add 50 $\mu$ L of spore suspension from spore concentration of  $8 \times 10^4$  spores/mL to make the final count of spores in the well as  $2 \times 10^4$  spores/mL.
17. Add 50 $\mu$ L of Nipagin to the triplicates.
18. For blank, add 100 $\mu$ L of PDB and 100 $\mu$ L of PBS.

Well1												
	1	2	3	4	5	6	7	8	9	10	11	12
A	Conc. 1	Conc. 2	Conc. 3	Conc. 4	Conc. 5	Conc. 6	Conc. 7	Conc. 8	Conc. 9	Conc. 10	Conc. 11	Conc. 12
B												
C												
D	Fungal spores with ethanol			Blank	Fungal spores with ethanol			Blank				
E	Only fungal spores				Only fungal spores							
F	Conc. 1	Conc. 2	Conc. 3	Conc. 4	Conc. 5	Conc. 6	Conc. 7	Conc. 8	Conc. 9	Conc. 10	Conc. 11	Conc. 12
G												
H												

Row A, B, C: 100 $\mu$ L PDB+50 $\mu$ L spore suspension of  $8 \times 10^4$  spores/mL (*Aspergillus versicolor*)+ 50 $\mu$ L Nipagin(Conc in  $\mu$ g/mL)

Row E, F, G: 100µL PDB+50µL spore suspension of  $8 \times 10^4$  spores/mL (*Rhizopus oryzae*)+ 50µL Nipagin(Conc in µg/mL)

D1, D2, D3: 100µL PDB+50µL spore suspension of  $8 \times 10^4$  spores/mL (*Aspergillus versicolor*)+ 50µL 100% ethanol

D7, D8, D9: 100µL PDB+50µL spore suspension of  $8 \times 10^4$  spores/mL (*Rhizopus oryzae*)+ 50µL 100% ethanol

E1, E2, E3: 100µL PDB+ 100µL spore suspension of  $4 \times 10^4$  spores/mL (*Aspergillus versicolor*)

E7, E8, E9: 100µL PDB+ 100µL spore suspension of  $4 \times 10^4$  spores/mL (*Rhizopus oryzae*)

19. Keep the plate in a shaker cum incubator at 28°C.

20. Take the OD measurements at 605nm in a plate reader for a fixed interval of time.. (For our purpose, we took 4 times a day at an interval of 6 hr for 4-5 days.)

21. Record the reading.